



LDMS User Support
Phone: 716-834-0900 x7311
E-mail: ldmshelp@fstrf.org
www.ldms.org/contact/

LDMS Training Team
Phone: 716-834-0900 x7466
E-mail: ldms.training@fstrf.org

The Reagent Logging module is used to track reagents used in the laboratory and uses a defined set of reagent codes. The module will record the Lot no., expiration date, manufacturer and start/end dates. After entered in the log, the lot number is available in the Edit Primary and Edit Aliquot menus to be assigned to a specimen containing the reagent.

Reagent Logging

Add Reagent

Reagent Name	Lot Number	Expiration Date	Manufacturer	Start Date	End Date	Active	In Use	Action
FBS	123456	10/Mar/2018	Fisher Scientific	06/Nov/2017	10/Jan/2018	<input checked="" type="checkbox"/>		Delete

Selected Reagent

Reagent Name: *

Lot Number: *

Expiration Date: *

Manufacturer: *

Start Date: *

End Date: *

Active:

1. From the LDMS menu bar, hover over **Administration** and select **Reagent Logging**.
2. **Click** the **Add Reagent** button at the top-right of the grid.
3. A blank like appears in the grid, highlight the line.
4. In the **Selected Reagent** menu, enter the required fields.
5. Load the entry into the grid by clicking, **Update Reagent**.
6. Click **Save**.

Selected Reagent Menu

Reagent Name: This field contains the three letter code for the reagent that is being tracked.

Manufacturer/ Lot Number: These fields are free text fields enter the manufacturer name and lot # from the container

Expiration Date: Manufacturer set date of expiration

Start/End Date: These fields record the date the reagent was first used and the date use was discontinued.

Active: New entries default to active.

To edit a previously entered lot. Select the entry in the grid. Adjust any of the fields in the **Selected Reagent** menu and click **Update Reagent**

Linking Reagent to Sample in Specimen Management

Aliquots for 0500-001ZLC00-000

Global Specimen ID	Status	Derivative Type	Sub Add/ Der Type	Specimen Condition	Available Volume	Other Specimen ID	Action
0500-001ZLC00-001	✓ [Icons]	PL2	N/A	SAT	1 ML		Edit [Dropdown]
0500-001ZLC00-002	✓ [Icons]	PL2	N/A	SAT	1 ML		Edit [Dropdown]
0500-001ZLC00-003	✓ [Icons]	PL2	N/A	SAT	1 ML		Edit [Dropdown]
0500-001ZLC00-004	✓ [Icons]	CEL	DMS	SAT	3000000 CEL		Edit [Dropdown]
0500-001ZLC00-005	✓ [Icons]	CEL	DMS	SAT	3000000 CEL		Edit [Dropdown]
0500-001ZLC00-006	✓ [Icons]	CEL	DMS	SAT	3000000 CEL		Edit [Dropdown]

Edit Aliquot

FRONTIER.0111111
 TEST
 6 Wk 01/Oct/2015
 0500-001ZLC00-000
 0500-001ZLC00-004

Global Specimen ID: 0500-001ZLC00-004

Specimen ID:

Other Specimen ID:

Enrollment: FRONTIER/TEST *

Derivative Type: CEL *

Sub Add/Der Type: Sub Add/Der Type

Sub Add/Der Reagent: *

Is Available:

After the lot is entered in the Reagent Logging module, the lot number is available in the **Edit Aliquot** menu to be applied to any sample using the sub additive.

- In **Specimen Management**, use the filters to navigate to the specimens page.
- In the **Aliquot grid**, locate the specimens utilizing the logged reagent
- Click the **Edit** button, the **Edit Aliquot** menu will open in a pop up window
- Select the appropriate lot number in the **Sub Add/Der Reagent** menu
- Click **Save**

Linking Additive to Reagent Lot

- In **Specimen Management**, click the **Edit** button for the Primary containing the logged Additive
- In the **Edit Primary Specimen** window select the lot number in the **Additive Reagent** menu

Additive:

Additive Type: EDT *

Additive Reagent:

Note: If known, the reagent lot number can be applied in Quick Add at the time of accessioning. In the **Aliquot grid**, locate entry with the Sub Additive/Derivative. Click the **Edit** button, select the lot number in the **Sub Add/Der Reagent** menu.